



D^r DOT

**BlueDot
Interpretation Software
for Windows**

**Dr DOT Software
allows quick and easy
Interpretation of**

**BlueDot
test results**

www.d-tek.be

Your Choice for a reliable serology in Autoimmunity and Gastroenterology






D^r DOT

BlueDot
Interpretation Software
for Windows

- dedicated to the interpretation of the whole range of BlueDot kits
- can be used with an office flatbed scanner
- quantification of the results allows an objective test interpretation
- contains a database that keeps record of all saved results (numeric values and scanned images)
- includes a multi-criteria research engine
- multilingual user friendly interface
- Results (values and scanned image) → can be sent by e-mail (Outlook or Outlook Express)
→ can be printed
→ can be exported in CSV format
- 30 days evaluation period free of charge
- lifelong access licence including free of charge updates



System requirements :

PC Windows/98/ME/2000/NT/XP with 800 Mhz processor, 64 Mb RAM and 100 Mb free space on hard disk. Flatbed A4 scanner compatible TWAIN with 300 dpi mode and 250 levels of gray, installed and configured appropriately. Recommended scanners: Epson Perfection series. Option: printer, internet connection.

