

**MATERIAL SAFETY DATA SHEET**



**Connectivitis<sup>10</sup> IgG**

**1. Product Identification**

|                      |  |
|----------------------|--|
| <b>Product Name:</b> | BlueDot Connectivitis <sup>10</sup> IgG  |
| <b>Product Code:</b> | CT10D-24   |
| <b>Intended Use:</b> | Immunodot kit for the detection of IgG antibodies to Nucleosomes, dsDNA, Histones, Sm, Sm/RNP, SSA/Ro60kD, TRIM21 (SSA/Ro52kD), SSB, Jo-1 and Scl-70 antigens in human serum |
| <b>Company Name:</b> | D-TEK s.a<br>Parc Initialis, rue René Descartes 19<br>BE-7000 Mons Belgium<br>Tel. : +32 65 841 888  |

**2. Composition, Information on ingredients**

• **Classification : Mixture**

|                      |                                     |
|----------------------|-------------------------------------|
| <b>Abbreviations</b> | <b>Substances</b>                   |
| AP                   | Alkaline Phosphatase                |
| BCIP                 | Bromo-Chloro-Indolyl-Phosphate      |
| BSA                  | Bovine Serum Albumin                |
| MIT                  | MethylIsoThiazolones (preservative) |
| NBT                  | NitroBlueTetrazolium                |
| TBS                  | Tris Buffer Saline                  |

|                     |   |
|---------------------|---|
| <b>Kit Contents</b> | <b>Composition</b>  |
| Sample Buffer       | TBS, Tween, BSA, MIT (15 ppm).  |
| Wash Buffer         | TBS, Tween, MIT (15 ppm).   |
| Conjugate           | Goat anti-human IgG/AP, MIT (15 ppm).   |
| Substrate           | NBT/BCIP, Sodium Azide (0.05 %)   |
| Membrane Strips     | Membrane-coated purified Nucleosomes (natural, bovine), dsDNA (natural, bovine), histones (natural, bovine, mixture of H1, H2a, H2b, H3 and H4), Sm (natural, bovine), Sm/RNP (natural, bovine), SSA/Ro60kD (human recombinant), TRIM 21 (SSA/Ro52kD, human, recombinant), SSB (recombinant, human), Jo-1 (recombinant, human) and Scl-70 (natural, rabbit) antigens. |

• **Hazardous Substances**

| Name         | CAS / EINECS         | Concentration in mixture | R / S phrase (at concentration used) | Symbol(s) (at concentrations used) |
|--------------|----------------------|--------------------------|--------------------------------------|------------------------------------|
| Sodium Azide | 26628-22-8/247-850-1 | < 0.1 %                  | -<br>(non hazardous)                 | -<br>(non hazardous)               |

Where an ingredient in the mixture is not listed as hazardous substance, this means that the substance has not been classified by the European Communities legislation on dangerous substances. This does not mean that it has no hazardous properties: please refer to further recommendations and precautions listed in this document for maximum safety when using this product.

• **Other information**

Bovine serum albumin (BSA) used in the preparation of reagents in this kit was derived from bovine blood collected at a USDA licensed establishment. The cattle received ante- and post-mortem health inspection at the abattoir by a U.S. Veterinary Services Inspector and they were apparently free from infectious and contagious diseases. All donor animals were sourced from the United States, a country in which Bovine Spongiform Encephalopathy is not known to exist.

**3. Health hazard data / Risk identification**

Liquid reagents contained in the kit may be irritating to eyes and skin

**Main exposure risks and symptoms**

|                           |   |
|---------------------------|---|
| <b>Contact with eyes:</b> | Irritation. Tears   |
| <b>Contact with skin:</b> | Irritation  |
| <b>Ingestion:</b>         | Although not considered particularly toxic or dangerous upon ingestion in conditions of normal use, avoid ingestion and contact with food |

**4. Emergency and first aid procedures**

|                           |  |
|---------------------------|--|
| <b>Contact with eyes:</b> | Immediately flush eyes thoroughly with water.  |
| <b>Contact with skin:</b> | Immediately wash skin with soap and large volumes of water.  |
| <b>Ingestion:</b>         | If swallowed, wash out mouth with water provided the person is conscious; seek medical advice (showing this document when possible). Never give anything by mouth to an unconscious person; never try to make an unconscious person vomit. |

**5. Fire hazard data**

|  |  |
|--|--|
| <b>Flammability:</b>                     | Liquid reagents contained in the kit are not combustible.<br>Combustion of cardboard inserts inside the kit and the outer cardboard box of the kit may produce intense heat.   |
| <b>Extinguishing Media:</b>              | Water, carbon dioxide, dry chemical powder or polymer foam.<br>Use extinguishing media appropriate to surrounding fire conditions.   |
| <b>Special Fire Fighting Procedures:</b> | For fires involving this material, do not enter any enclosed or confined fire space without proper protective equipment. This may include self-contained breathing apparatus to protect against the hazardous effects of the normal products of combustion or oxygen deficiency. |

## MATERIAL SAFETY DATA SHEET

### 6. Accidental release measures

|                                      |  |
|--------------------------------------|--|
| <b>Precautionary measures:</b>       | Wear appropriate protective clothing (refer to & 8)  |
| <b>Environmental precautions:</b>    | Avoid flushing away in drainpipes.   |
| <b>Released or spilled material:</b> | Sweep up and collect in appropriate container for waste disposal (refer to & 13); clean the floor and all other contaminated objects with water. |

### 7. Handling and storage

|                                |  |
|--------------------------------|--|
| <b>Handling:</b>               | Wear appropriate protective clothing (refer to & 8) ; eliminate rapidly from eyes, skin and clothes ; wash hands and any other exposed zone with water and mild soap before eating, drinking, smoking and leaving workplace. |
| <b>Storage:</b>                | Store in a cool dry place  |
| <b>Additional information:</b> | refer to & 10  |

### 8. Exposure controls and personal protection

|                          |                                       |
|--------------------------|---------------------------------------|
| <b>Hands protection:</b> | Laboratory gloves                     |
| <b>Skin protection:</b>  | Laboratory coat                       |
| <b>Eyes protection:</b>  | Safety goggles                        |
| <b>Ingestion:</b>        | Do not eat, drink or smoke during use |

### 9. Physical and chemical properties

|                |   |
|----------------|---|
| <b>Form:</b>   | Liquid reagents   |
| <b>Colour:</b> | colourless (Wash), yellow (Diluent), red (Conjugate), pale yellow (Substrate) |
| <b>Smell:</b>  | negligible  |

### 10. Stability and reactivity data

|  |  |
|--|--|
| <b>Dangerous reactions:</b>                            | Particular dangerous reactions not known   |
| <b>Dangerous properties:</b>                           | Particular dangerous properties not known  |
| <b>Hazardous combustion or decomposition products:</b> | Nature of decomposition products from liquid reagents: not known.<br>Combustion of cardboard inserts inside the kit and of the outer cardboard box of the kit does <u>not</u> liberate toxic gas (only carbon dioxide and water vapour). |

### 11. Toxicological information

Data not yet available

### 12. Ecological information

Data not yet available

### 13. Disposal considerations

Contact a licensed professional waste disposal service to dispose of this material. Observe all federal, state and local environmental regulations. Emptied bottles and vials may retain product residues : handle as if they were full.

### 14. Transport information

Non classified

### 15. Regulatory information

From the European Communities legislation on dangerous substances up to and including EC Council Directive 96/54/EC, this product is not hazardous and does not contain hazardous ingredients at concentrations used.

### Other information

Contents and format of this Safety Data Sheet comply with the EC Council Directive 93/112/EC.

The above information is believed to be correct but does not purport to be all inclusive and shall be used only as a guide. D-TEK s.a. shall not be held liable for any damage resulting from handling or from contact with the above product.